







Easy to take apart & reassemble

Don't think - do it! Less than a minute is all you need to take your *Invacare® Colibri™* apart ready for transportation.

Thanks to the ingenious and unique *Invacare® LiteLock™* system the scooter dismantles easily and fits into any car boot. No tools are needed, and the click-together mechanism makes everything quick and easy to take apart and reassemble, staying securely in position.

So drive it anywhere your fancy takes you. The convenient, spacious front basket can store everything you pick up along the way.

Moreover, the scuff-free wheels won't leave any unsightly black marks behind you.









The *Invacare*[®] *Colibri*[™] Mix n' Match your style, with ease





Mix n' Match your style, with ease

Why should a scooter be the same colour every day? Why can't you coordinate it to your clothes or to simply match your mood?

The *Invacare® Colibri™* is designed with the latest trends in mind, so you can pick the colours you like best and brighten up the day with your preferred shades!

Colour panels can easily be removed and clipped back on without the need for any tools. The panels' clip-on mechanism is quick and easy, guaranteeing everything fits together perfectly and stays fixed for as long as you wish - until you fancy a change.











The *Invacare*[®] *Colibri*[™] Safe, manoeuvrable & comfortable





Safe, manoeuvrable & comfortable

The *Invacare® LiteLock™* system avoids any mistakes during the dismantling and assembly of your scooter. Moreover it easily brings the scooter's parts together and ensures they securely stay in position, thus delivering a safe driving experience.

Extra built-in reflectors have been added to further increase your safety.

The *Invacare®* **Colibri**TM is designed to be one of the most portable micro scooters available. With its small footprint the **Colibri**TM offers greater manoeuvrability in tight, limited spaces.

An ergonomic handlebar offers more hand position options, helping to prevent fatigue in the hands and wrists. The comfortable padded seat with elegant colour stitching also ensures a smooth relaxing drive at all times.













Personalise your Invacare® Colibri™

Seat suspension for more comfort

For increased comfort during outdoor driving add the optional seat suspension kit. The robust seat suspension makes a big difference over uneven terrains.



Choose from 12 or 18 Ah battery packs

For those who like to travel further with peace of mind, there is the option of a larger battery pack. Plus, charging the battery is easy — either on the scooter itself, or detached, in the garage or your home.



Storage cover

If you store your scooter outside you can protect it from the elements by adding a storage cover when not in use. This will keep it dry, protecting the workings of the scooter.



Crutch & cane holder

For those who require a crutch or a cane, keep it securely stored until you reach your destination.



Wheels

The *Invacare*® *Colibri*TM comes with either smaller wheels allowing you to go in confined spaces or with a bigger wheel set for those who prefer to drive outdoor in more uneven terrain, on cobbled streets.



Personal key

For extra safety you can secure your scooter by adding a personalised key. As the cylinder initially comes with just two keys, we recommend that you make some copies for added security and peace of mind.



Rear basket


The rear basket is perfect for shopping, and can be installed without any tools. Just attach it when needed and store away once the shopping's done!





The Invacare® Colibri™
The Colourful route to effortless, safe driving





Technical data



465 mm



400 mm



370 - 420 mm



385 mm



505 mm



1010 mm



44 kg (12 Ah batteries) 49 kg (18 Ah batteries)



136 kg



2 x 12 Ah 2 X 18 Ah



2200 mm



45 mm



6° / 10.5%



200 W



Up to 16 km



8 km/h

Dealer's stamp

